

Using uDirect 3 with a mobile phone

If you have a mobile phone equipped with Bluetooth functionality (e.g. iPhone, Android phone), uDirect 3 may be used to stream phone conversations wirelessly to your hearing aids allowing you to hear the callers voice directly in your hearing aids. uDirect 3 picks up your voice through its own microphone and streams it to your phone.



During phone calls, the access (A) button provides a range of features that help you manage your calls.

To benefit from the pairing of the uDirect 3 and your phone, we recommend wearing the uDirect 3 around the neck. During calls, your voice can only be picked up if the uDirect 3 is not covered by clothes.






Pairing uDirect 3 with your mobile phone

1. Enable Bluetooth on your phone settings menu.
2. Ensure that the uDirect 3 is turned ON and the neckloop is plugged in and worn correctly. Press the connect (↔) and volume up (+) buttons at the same time for two seconds. The audio indicator will blink blue to indicate pairing mode on the uDirect 3.
3. Enable Bluetooth pairing mode on your phone so that a new device can be added. See manufacturer instructions for details. If prompted, enter “0000” for pairing code.
4. If prompted, accept request to grant uDirect 3 access to your phone book. This is used for voice calling and caller ID features.
5. After successful pairing, the blue audio indicator on the uDirect 3 will stop blinking.



Once pairing is done, your phone and uDirect 3 should stream audio automatically.

Notes:

-  Please contact your hearing healthcare professional if you can not complete the pairing procedure with the steps described.
-  To receive phone calls your phone must support either the Bluetooth “Head Set” (HSP) or “Hands Free” (HFP) profile.
-  To access caller ID features, your phone must support the Phone Book Access Profile (PBAP)