



Element™ is a series of high-quality, brilliant sounding hearing instruments featuring premium innovations, a range of automatic functionality, and a personalized fit in the advanced through essential categories.

element at-a-glance

Performance Profile

	e ¹⁶	e ⁸	e ⁴
Range of loss	Mild to severe	Mild to severe	Mild to severe
Channels	16	8	4
Bands	16	8	8
Processing type	WDRC Linear Limiting	WDRC Linear Limiting	WDRC Linear Limiting

Features

Automatic programs	AutoPro2™	AutoMic™	
Manual listening programs	up to 3	up to 3	up to 3
Ability to disable listening programs	•	•	•
Directional microphone system	Adaptive directional, directional, omni	Adaptive directional, directional, omni	Directional, omni
AntiShock™	•	•	•
Speech enhancement	Speech Enhancement LD	•	
Noise reduction	•	•	•
Wind noise manager	•	•	•
Feedback management	Phase canceller	Phase canceller	Phase canceller
Automatic telephone solution	Easy-t (Custom and BTE)	Easy-t (Custom and BTE)	Easy-t (Custom and BTE)
Data logging	•	•	•
Slim tube/open fitting	Moda and 13 no VC	Moda and 13 no VC	Moda and 13 no VC
Quiet mode expansion	•	•	•
Manual volume control (may be disabled)	•	•	•
Start up delay	•	•	•
Ideal volume indicator	•	•	•
Program beep indicator	•	•	•
Low battery warning	•	•	•
Optional telecoil	•	•	•
MLx compatible	BTE (except Moda)	BTE (except Moda)	BTE (except Moda)

this is your element



www.unitronhearing.com