



You can have it all  
with a little Moxi.™



moxi moxi moxi moxi moxi moxi moxi moxi moxi moxi  
moxi moxi moxi moxi moxi moxi moxi moxi moxi moxi

So much to enjoy. So little to see.



Discreet, attractive,  
lightweight style

Digital technology smart  
enough to adjust to your  
listening needs

Same-day fitting lets  
you hear what you  
want instantly

Moxi style is available at four  
levels of technology designed  
for your lifestyle, listening  
preferences and budget

# brilliant

Moxi™: designed to deliver brilliant sound.

## Speaker unit

Thin tube that discreetly connects the hearing device to the speaker

## Speaker

Rests in your ear canal enhancing sound

## Microphone

Automatically adjusts to focus on what you want to hear

## Push button

Gives you control over sound performance in your environment

## Battery door

Offers instant and easy access to change your battery

## Control when you need it:

Unique to such a small, discreet hearing device, the push button helps you control your volume or programs in diverse listening situations.



Look  
closely.

Look

Moxi's innovative, sleek design is 37% smaller than the average Behind-The-Ear (BTE) device. Look closely – you'll be sure to love what you see, and what you can't see.



**What is a BTE?**

A BTE (Behind-The-Ear) is a hearing device that fits snugly behind your ear.

# closely



## More than meets the eye.

Strikingly discreet, with a virtually invisible tube, no one will know Moxi's there. Yet, it's so attractive you might be tempted to show it off.

## Go ahead. Pick a color.

It's your choice: make a fashion statement with Moxi's signature black hearing device or choose one that matches with your skin or hair color.



Beige



Brown



Light Grey



Dark Grey



Taupe



Black



Pearl



# big **Small size, big comfort.** comfort



## **Created for comfort.**

Moxi offers a lightweight, open-fit design that leaves your ear more open so you can hear more naturally. Delivering premium comfort and brilliant sound that makes all-day wear more enjoyable.

## **Fit for your lifestyle.**

Your life is full of activity. It's how you choose to live. When you select Moxi, you're instantly fit with your new digital hearing device. That's right, the very same day.



# natural sounding experience.

A natural

sounding experience.



**So much technology.  
Such a little Moxi.**

Moxi's cutting-edge canal receiver technology (CRT) uses a thin tube to connect with a speaker in your ear canal. Unlike a conventional BTE where the speaker is inside your hearing device, with CRT the speaker is actually in your ear producing a more brilliant sound experience.

**Moxi: sound the way it  
was meant to be heard.**



speaker

**Conventional  
approach**



speaker

**Canal Receiver  
Technology  
approach**

# Don't miss a single moment.

## Hear what matters most.

Music plays softly, serving spoons clatter in dishes, while you chat with dinner guests. You nod in agreement, but only catch part of what's being said. Imagine a scenario where

you easily understand conversation – feel comfortable in diverse listening situations – and can connect with those who matter most. That's your life with Moxi.

**It's not about hearing.  
It's about understanding  
everyday conversation.**



# a moment



## Live life with Moxi:

- ▶ improved speech understanding
- ▶ ability to hear crisp, natural sound
- ▶ enjoying conversation in a variety of listening environments
- ▶ participating in more social and leisure activities
- ▶ reduced listening fatigue
- ▶ renewed confidence

# Reconnect

Reconnect with your  
**surroundings.**

## **Automatic performance**

Moxi's automatic performance invites you to move seamlessly between situations – from home to work, intimate conversation to family dinner. As sound changes, Moxi adapts to your needs providing the best listening experience.

## **Directional microphones**

Moxi's directional microphones produce clear, intelligible speech by emphasizing the sound in front of you. When you change listening environments, Moxi automatically selects the most effective microphone setting, giving you the confidence to conduct a conversation anywhere.



# Everyday

Everyday sound and noise  
**management.**

## **Speech enhancement**


Moxi increases the amplification of speech and improves your understanding of every nuance of a conversation.

## **Noise reduction**

Moxi focuses on speech without over-amplifying background noise like the sound of a dishwasher. It gives you the support to participate in any situation.

## **AntiShock™**

A cool little feature that's smart enough to know when not to amplify noises, like a car door slamming, that could be uncomfortable. Moxi delivers instant comfort in a variety of listening situations.



# So many cool

So many cool things to know.

## All in a little Moxi.™

**Open the lines of communication.**

Moxi automatically switches to its telephone solution making it even more enjoyable to chat with an old friend.



# things

## **The perfect pair.**

An optimal hearing solution starts with two hearing devices. When sound is equally amplified in both your ears, Moxi can support you anytime and anywhere.



## **Sound the way you want it.**

Hear what you want at a noisy restaurant or gathering. Moxi gives you the confidence to participate in conversations no matter where you are.

# Get some **Get some** **sound advice.**

## **Make the right choice together.**

Every hearing loss is unique. Your hearing healthcare professional is the most qualified person to help select the right Moxi for you. Together, you'll address your listening challenges and preferences, ensuring maximum performance from your Moxi.

## **A product for every lifestyle.**

Whether you have an active life or a quiet one, we've got the Moxi for your lifestyle, listening preferences and budget. Select from these four unique products styled with Moxi: Indigo™, Element™ 16, Element 8 and Element 4.

**4 products. 1 is right for you.**



# sound advice



	indigo moxi	e <sup>16</sup> moxi	e <sup>8</sup> moxi	e <sup>4</sup> moxi
<b>My listening situations</b>	Busy, active life, many different listening environments	Active life, multiple environments	Social time spent at home and in other diverse environments (e.g. restaurants)	Social time spent in fairly quiet, less demanding situations (e.g. watching TV, playing cards)
<b>My priorities</b>				
Reaction to changing environments				
Premium innovative technologies				
Comfort and convenience technologies				
<b>My budget</b>				

moxi moxi moxi moxi moxi moxi moxi moxi  
moxi moxi moxi moxi moxi moxi moxi moxi



### **About Unitron Hearing**

Unitron Hearing creates high-value, easy-to use hearing solutions that make everyday listening easier for people with hearing loss. The Company works with hearing healthcare professionals around the world to deliver innovative hearing devices that meet the needs of every hearing loss and lifestyle.

