


Unparalleled fitting flexibility.

Only Moxi offers both a standard xReceiver and a power xReceiver within one CRT product family, to satisfy clients with hearing losses ranging from mild to severe. A wide array of coupling and venting options ensures optimal performance for every client. Moxi brings instant fit capability or custom options to both standard and power fittings.

Standard xReceiver



Open dome:
up to 65 dB HL

Choose for:

- ▶ Instant fit
- ▶ Open fit

Considerations:

- ▶ Occlusion free



Closed dome:
up to 70 dB HL

Choose for:

- ▶ Instant fit
- ▶ Semi-occluded fit

Considerations:

- ▶ Extended fitting range over open dome
- ▶ Some occlusion



Power dome:
up to 75 dB HL

Choose for:

- ▶ Instant fit
- ▶ Occluded fit

Considerations:

- ▶ Extended low-frequency response
- ▶ Occlusion



Sleeve mold:
up to 75 dB HL


Choose for:

- ▶ Custom fit
- ▶ Choice of venting

Considerations:

- ▶ Comfort
- ▶ Possibly occlusion, depending on venting
- ▶ Impression required

Power xReceiver



Power dome:
up to 95 dB HL

Choose for:

- ▶ Instant fit



Power shell:
up to 95 dB HL

Choose for:

- ▶ Custom fit
- ▶ Comfort
- ▶ Choice of venting

Considerations:

- ▶ Impression required