



Remote adjust

Remote adjust gives you the ability to make frequency response and adaptive feature changes in TrueFit fitting software and remotely send the changes along with a personalized message to your clients through the Remote Plus app*.



A Sonova brand

Meet expectations without meeting in-person

Barriers such as distance, accessibility, schedules and health concerns no longer stand in the way of providing professional care that's practical for both you and your clients. Now you can provide a comprehensive and convenient distance support solution that doesn't require them to leave home, or you to use up an appointment slot.

A continuous care feature that's not confined to the in-clinic appointment:

- Changes are applied remotely without the client physically present in the clinic
- Adjustments can be made to the following parameters in any automatic, manual or optional app program: gain offset, microphone mode, dynamic noise reduction, speech enhancement, noise reduction, AntiShock 2, phase canceller, and wind control
- No need to virtually connect or sync with the client – you submit the changes, and they can apply them at any time



You experience ...

- The ability to service clients that need minor adjustments quickly and efficiently
- An opportunity to free up appointment slots



Clients experience ...

- Convenience of receiving program adjustments without leaving home
- Reassurance they can revert back to previous settings if changes don't meet their needs

Watch the videos today by scanning the QR code!

We have created a 5 minute self-paced online learning session to provide an overview of Remote Adjust. To watch this webinar, you will need to login with your myUnitron account information. You will then see the icon to watch this webinar series. If you are not sure about your account information, please contact Jerome Wheatcroft at Jerome.Wheatcroft@unitron.com or at 1.800.265.8255.



*For the use of the Remote Plus app, Unitron hearing aids with Bluetooth® connectivity are required. The Remote Plus app works on Apple smartphones with iOS 12 or newer and Android smartphones with version 7 or newer. App Store badge is a trademark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC.