

Unitron's Innovative Product Family

LET'S HAVE A TOAST TO
*hearing what
matters most*

Welcome to the
Vivante.
experience

Introducing Vivante, Unitron's newest platform, built to bring your clients the sounds of the good life right to their ears. HyperFocus for face-to-face conversations in loud noise, a 360 conversation in car manual program that focuses on speech regardless of where it's coming from, new RIC and charging case designs, all delivering a personalized listening experience.



Receiver-in-canal (RIC) Mild to severe



Moxi V-R
(Li-ion rechargeable)



Moxi V-RT
(Li-ion rechargeable with telecoil)

Battery size

Li-ion

Available in

9 7 5 3

Features


- New design
- Rechargeable with new charger
- Tap control
- Made for all direct connectivity
- Connect to two different Bluetooth® enabled devices at once

- Phone and media connectivity to both ears
- Multi-function buttons
- Choice of standard, moderate, power and ultra power receivers
- Telecoil (Moxi V-RT)
- 2.4 GHz wireless
- Plasma coating
- IP 68 Ingress Protection rating


Processing type

- WDRC
- Linear limiting


Fitting range




Standard Receiver (S)
(111/47)



Moderate Receiver (M)
(114/51)



Power Receiver (P)
(122/59)



Ultra Power Receiver (UP)
(132/71)

— Open dome/cap dome
- - - Vented dome
... Power dome/SlimTip

Moxi

Receiver-in-canal (RIC) Mild to severe

Stride

Behind-the-ear (BTE) Mild to severe

Powered by **Blu**



Moxi B-R
(Li-ion rechargeable)



Moxi B-312



Moxi B-RT
(Li-ion rechargeable with telecoil)

Battery size

Li-ion (Moxi B-R, Moxi B-RT)
312 (Moxi B-312)

Available in

9 7 5 3



Features

- Tap control (rechargeable only)
- Made for all direct connectivity
- Connect to two different Bluetooth enabled devices at once
- Phone and media connectivity to both ears
- Rechargeable, Li-ion (Moxi B-R, Moxi B-RT)
- Multi-function buttons
- Choice of standard, moderate, power and ultra power receivers
- Telecoil (Moxi B-RT only)
- 2.4 GHz Wireless
- Plasma coating
- IP 68 Ingress Protection rating

Processing type



- WDRC
- Linear limiting

Fitting range



Standard Receiver (S)
(111/47)

Moderate Receiver (M)
(114/51)



Power Receiver (P)
(122/59)

Ultra Power Receiver (UP)
(132/71)

— Open dome/cap dome
- - - Vented dome
... Power dome/SlimTip



Stride™ B-PR
(Li-ion rechargeable)



Stride™ B-UP

Battery size

Li-ion (Stride B-PR)
675 (Stride B-UP)

Available in

9 7 5 3


Features

- Made for all direct connectivity
- Phone and media connectivity to both ears
- Rechargeable, Li-ion (Stride B-PR)
- Multi-function buttons
- Earhook or slim tube
- 2.4 GHz Wireless
- Plasma coating
- IP 68 Ingress Protection rating
- Telecoil (Stride B-UP only)

Processing type


- WDRC
- Linear limiting

Fitting range



..... Slim tube (69/129)

Stride B-PR
(130/68)



Stride B-UP
(141/84)

Insera

Full shell/half shell/canal/mini-canal/CIC/IIC

Battery size

10A
312

Available in

9 7 5 3

Features

- Made for all direct connectivity*
- Connect to two different Bluetooth enabled devices at once*
- Virtually invisible in the ear (IIC)
- Choice of 3 power levels (Omni)
- Choice of 4 power levels (B-312)
- IntelliVent
- Easy-t

- Telecoil
- Wireless*
- Directional option*
- Push button battery door (option for IIC)**
- Remote functionality†
- Plasma coated door, guards, push button and VC (IIC door and push button only)

Processing type

- WDRC
- Linear limiting

*Only for Insera B-312

**CIC and mini canal only

†312 and IIC only



Insera™
B-10 NW O




Insera
B-312 NW O




Insera
B-312


Faceplate colors




Fitting range




Moderate (M)
(109/40)



Super Power (SP)
(119/60)



Power (P)
(115/50)



Ultra power (UP)
(127/70) 312 wireless only

Stride

Powered by **Discover Next**

Behind-the-ear (BTE) Mild to severe



Stride™ M

Battery size

312 (Stride M)

Available in

9 7 5 3

Features


- Made for all direct connectivity
- Phone and media connectivity to both ears

- Multi-function buttons
- Earhook or slim tube
- Telecoil
- 2.4 GHz Wireless
- Plasma coating
- IP 68 Ingress Protection rating

Processing type

- WDRC
- Linear limiting

Fitting range



Stride M
(128/63)

Insera

Full shell/half shell/canal/mini-canal/CIC/IIC

Battery size

10A
312

Available in

9 7 5 3

Features

- Virtually invisible in the ear (IIC)
- Choice of 3 power levels (Omni)
- Choice of 4 power levels (312 W Dir)
- IntelliVent
- Easy-t
- Telecoil

- 2.4 GHz Wireless *
- Directional option*
- Push button battery door (option for IIC)**
- Remote functionality†
- Plasma coated door, guards, push button and VC (IIC door and push button only)

Processing type

- WDRC
- Linear limiting


*312 W Dir only

**CIC and mini canal only


†312 Wir Dir and IIC only



Faceplate colors



Fitting range



Earmolds

Earmold styles

 Skeleton (Acrylic, silicone)	 Full Shell (Acrylic, silicone)	 Half Shell (Acrylic, silicone)	 Semi-Skeleton (Acrylic, silicone)	 Canal (Acrylic, silicone)	 Canal Lock (Acrylic, silicone)
 CROS (Acrylic, silicone)	 Hollow Sleeve (Acrylic)	 Solid Sleeve (Acrylic, silicone)	 cShell (Acrylic)	 RIC Secure Fit (RSF) (Silicone)	Not shown: uncarved full shell

Translucent (Acrylic & Silicone) Earmolds

 Trans. clear (21)	 Trans. pink (T)	 Trans. brown (N)
---	--	---

Solid Silicone Colors

 White (19)	 Yellow (20)	 Orange (11)	 Red (10)	 Green (17)	 Blue (07)	 Purple (08)	 Black (06)	 Jungle Green (2B)	 Flamingo Pink (5B)
--	--	--	---	---	--	---	---	--	---

Metallic Silicone Colors

 Galactic Blue (1B)	 Star Dust Silver (3B)	 Magic Pink (4B)
--	--	--

Vivante features at-a-glance

Integra OS		9	7	5	3
Environmental classification	Total listening environments	8	6	4	2
	Conversation in noise	•			
	Conversation in a crowd	•			
	Conversation in a small group	•	•		
	Music	•	•		
	Noise	•	•	•	
	Quiet	•	•	•	
	Conversation in quiet	•	•	•	•
	Conversation in noise	•	•	•	•
	Total streaming environments	2	2	2	2
	Media speech	•	•	•	•
	Media music	•	•	•	•
Sound optimization and performance	HyperFocus	•	•		
	Conversation in loud noise (manual program)	•	•		
	360 conversation in car (manual program)	•			
	AutoFocus 360	•			
	AutoFocus		•		
	Dynamic noise reduction	•	•		
	Soft speech lift	•	•	•	
	Pinna Effect 2	•	•	•	
	Pinna Effect	•	•	•	•
	Speech enhancement	•	•	•	•
	Noise reduction	•	•	•	•
Hardware features	Tap control	•	•	•	•
	Made for All direct connectivity	•	•	•	•
	Rechargeable	•	•	•	•
	Telecoil (RT only)	•	•	•	•
Fine-tuning channels					
		20	20	16	12

Available in all technology levels

Sound stabilization

- › AntiShock 2
- › Wind control
- › Feedback manager

Ease and convenience

- › Remote Plus app
- › TV Connector
- › Remote Control
- › PartnerMic
- › RogerDirect
- › DuoLink
- › Binaural Phone

Microphone options

- › Adaptive directional
- › Fixed directional
- › Fixed wide directional
- › Omni directional

Fitting

- › First fit workflow
- › Automatic Adaptation Manager
- › Frequency compression 2
- › Tinnitus masker
- › Manual programs
- › IntelliVent

Experience Innovations

Performance tracking

- › FLEX:TRIAL™
- › Log It All
- › Ratings


Continuous care

- › Coach
- › Remote adjust
- › FLEX:UPGRADE™

Personalization

- › Equalizer
- › Clarity and comfort buttons
- › Optional app programs

Vivante Moxi RIC colors



Blu features at-a-glance

Integra OS		9	7	5	3
Environmental classification	Total listening environments	7	6	4	2
	Noise	•			
	Conversation in noise	•	•		
	Music	•	•		
	Conversation in a small group	•	•	•	
	Quiet	•	•	•	
	Conversation in a crowd	•	•	•	•
	Conversation in quiet	•	•	•	•
	Total streaming environments	2	2	2	2
	Streaming speech	•	•	•	•
	Streaming music	•	•	•	•
Sound optimization and performance	AutoFocus 360	•			
	AutoFocus		•		
	Dynamic Noise Reduction	•	•		
	Soft speech lift	•	•	•	
	Spatial Awareness	•	•	•	
	Pinna Effect	•	•	•	•
	Speech enhancement	•	•	•	•
	Noise reduction	•	•	•	•
	Microphone features	•	•	•	•
Hardware features	Tap control (PR, R and RT only)	•	•	•	•
	Made for All direct connectivity	•	•	•	•
	Rechargeable (PR, R and RT only)	•	•	•	•
	Telecoil (UP and RT only)	•	•	•	•
Fine-tuning channels					
		20	20	16	12

Available in all technology levels

Sound stabilization

- › AntiShock 2
- › Wind control
- › Feedback manager
- › Natural Sound Balance

Ease and convenience

- › Remote adjust
- › Remote Plus app
- › TV Connector
- › Remote Control
- › PartnerMic
- › RogerDirect
- › DuoLink
- › Binaural Phone
- › Tap control

Microphone options

- › Adaptive directional
- › Fixed directional
- › Fixed wide directional
- › Omni directional

Fitting

- › First fit workflow
- › Automatic Adaptation Manager
- › Frequency compression 2
- › Tinnitus masker
- › Manual programs
- › IntelliVent
- › MyMusic
- › Ultra Power pre-sets (UP only)

Experience Innovations

Performance tracking

- › FLEX:TRIAL™
- › Log It All
- › Ratings
- › Usage


Continuous care

- › Coach
- › Remote adjust
- › FLEX:UPGRADE™

Personalization

- › Equalizer
- › Clarity and comfort buttons
- › Optional app programs

Blu Moxi RIC and Stride BTE colors



Discover Next features at-a-glance

SoundCore		9	7	5	3
Environmental classification	SoundNav 3.0	•	•	•	•
	Total listening environments	7	6	4	2
	Conversation in a crowd	•			
	Conversation in a small group	•	•		
	Music	•	•		
	Quiet	•	•	•	
	Noise	•	•	•	
	Conversation in quiet	•	•	•	•
	Conversation in noise	•	•	•	•
	Total streaming environments	2	2	2	2
	MediaNav speech	•	•	•	•
	MediaNav music	•	•	•	•
Sound optimization	Sound Conductor	•	•	•	•
	Speech enhancement	•	•	•	•
	Soft speech lift	•	•	•	
	Noise reduction	•	•	•	•
	Microphone features	•	•	•	•
Localization	Spatial Awareness	•	•	•	
	Personalized	•			
	Pinna Effect	•	•	•	•
Performance in challenging environments	SpeechPro	•			
	Speech Locator	•	•		
	Speech Focus	•			
	Spatial Speech	•			
	Spectral Speech	•			
	SpeechZone 2		•		
Fine-tuning channels					
		20	16	14	12

Available in all technology level

Ease and convenience

- › Made for all direct connectivity
- › Phone and media connectivity to both ears (in some models)
- › Rechargeable*
- › TV Connector
- › DuoLink
- › Binaural Phone

Experience Innovations

- › FLEX:TRIAL™
- › FLEX:UPGRADE™
- › Insights
 - Remote adjust
 - Equalizer
 - Coach
 - Log It All
 - Ratings (including proactive alerts)
 - Overall satisfaction
 - Notifications
 - Usage
 - Data logging


Fitting

- › First fit workflow
- › Automatic Adaptation Manager
- › Frequency compression 2
- › Tinnitus masker
- › MyMusic
- › Manual programs
- › Streaming programs
- › NAL-NL2/NL1 and DSLv5
- › IntelliVent

Sound stabilization

- › AntiShock 2
- › Wind Control
- › Feedback manager
- › Natural sound balance

Discover Next Moxi RIC and Stride BTE colors



Accessories

Vivante, Blu and Discover Next products



Remote Control



PartnerMic



TV Connector



Unitron IIC remote



Li-ion RIC charger



Li-ion RIC Power pack



Li-ion BTE charger



Li-ion RIC Vivante charger



Remote Plus app*

roger

RogerDirect™


For more help in challenging situations like meetings or conversations in noise or over distance.

Compatibility may vary by product style.


*The Remote Plus app work on Apple smartphones with iOS 10.2 or newer and Android smartphone with version 6 or newer.

myUnitron


Online ordering, all your account activity—and more—right at your fingertips.



myStore – order products, accessories, earmolds, in-office materials and more



myServices – check order status, track shipments, warranty lookup, invoice copies, pay bills, software downloads and more



myLearning – online training and CEUs, call us today to get activated

24 hours a day. 7 days a week.

Register today at myUnitron.com or give us a call at 800.888.8882



Get ready to
Love the experience

At Unitron, we empower people with life-enhancing hearing experiences that fit seamlessly into their world. Our sound performance technology, experience innovations, and intuitive design work perfectly together for unmatched personalization and optimization. Because everyone deserves to **Love the experience.**[™]

800.888.8882
unitron.com/us

028-6248-03/03-23/hc © 2023 Unitron. All rights reserved.

sonova
HEAR THE WORLD