

# Unitron Custom Earpiece 4.0 Order Form

unitron™

Items on this order form are only compatible with Discover, Discover Next and Blu platform devices.

## Step 1: Order details

Account Number : \_\_\_\_\_

Hearing Clinic: \_\_\_\_\_

Delivery Address: \_\_\_\_\_

Contact Name: \_\_\_\_\_ Phone Number: \_\_\_\_\_

Email Address: \_\_\_\_\_

Date Ordered: \_\_\_\_\_ Date Required: \_\_\_\_\_

VIP 24 hour service (\$57.75 incl GST)

## Step 2: Patient information

Name: \_\_\_\_\_

Audiogram: (Required for Intellivent option)

kHz	0.25	0.5	1	2	4
LEFT AC (dB HL)					
RIGHT AC (dB HL)					

## Step 3: Product, style and options



SlimTip 4.0  
L R



cShell 4.0  
L R

Please note if a different configuration is required for each ear, please use two separate forms.

	Connection type				Shell colour				Vent style		Shell attribute	
	S Receiver	M Receiver	P Receiver	UP Receiver*	Clear (21)	Rose Tint (T)	Beige/Pink (26)	Tan (14)	Cavity	Regular	Hollow	Solid
SlimTip Acrylic				-	<input type="checkbox"/> S				<input type="checkbox"/> S		<input type="checkbox"/> S	-
SlimTip Silicone				-	<input type="checkbox"/> S		-	-	-	<input type="checkbox"/> S	-	<input type="checkbox"/> S
cShell Acrylic							<input type="checkbox"/> S		-	<input type="checkbox"/> S	-	-

\*HF3 filter recommended (recommended wax system for cShell xUP)  
S = Standard, will default if no other selection ticked

### SlimTip Acrylic options

Wax system ☐ Cerustop (CS) Ext. receiver tube + wax spring (WE) None

### SlimTip Silicone options

Wax system ☐ None

### Canal length

Short (RS) Medium (RM) - Standard Long (RL)

### cShell Acrylic options (HF3 recommended for UP receivers)

Wax system Cerustop (CS) HF 4 (H4) HF 3 (H3)\* Ext. receiver tube + wax spring (WE)  
Faceplate colour ☐ Pink (26) Cocoa (22) Brown (28) Transparent (13) Tan (14)

### Generic options (all earpieces)

Tube length (0-3) Left: Right:  
Vent size Intellivent - let Unitron choose. Occluded 1mm 1.1-1.6mm 1.7-2.2mm 2.3-3.0mm Open  
Removal line Transparent (21) Removal line extension + 5mm  
Locks Canal Skeleton

## Step 4: Preferences

If necessary, may we change the following:

Special instructions:

Step down power level - if audiogram permits (maintain 6dB Gain Reserve) Yes No

## Additional information

### Intellivent

We recommend to select Intellivent. We will set the correct vent diameter and style consider audiogram / shell geometry / feedback canceler and required low frequency gain. The client's audiogram is mandatory for the calculation of the Intellivent.

### Canal length

**Full length of impression** – The hearing instrument or earmould can be made the entire length of the impression. This is suggested for patients with moderate to profound hearing losses.


**Cut to mark** – This canal length is customised for each patient based on the mark you provide. The mark helps Unitron know exactly where to cut the canal for the individual patient's needs. This is suggested if a specific canal length is required or if the patient has soreness or a growth in the ear canal or a surgical ear.

**Extra-long** – The canal is cut 2 to 4 millimetres beyond the second bend. An extra-long canal is suggested for moderate to profound hearing losses.

**Long** – The canal is cut beyond the second bend. A long canal is suggested for moderate to profound hearing losses or for patients with occlusion.


**Medium** – The canal is cut between the first and second bends. A medium canal length is suggested for mild to severe hearing losses.

**Short** – The canal is cut just beyond the first bend. A short canal is suggested for mild to moderate hearing losses, narrow ear canals and/or patients with occlusion.



### Tube length

Measure the tube length to select the correct SlimTube. Measure both ears for binaural fitting.



### Hearing instrument compatibility

The SDS 4.0 earpiece portfolio is compatible with Discover or Discover Next hearing aids. It is not backward compatible to earlier generation hearing aids.

### Preferences

Let us know if you have specific requests. Select if hearing aid should be delivered together with the ear piece.