

# Accessory compatibility guide

Effective date: September 2018

Product	Tempus platform	North platform	Era platform
<b>For patients:</b>			
<a href="#">Rechargeable kit<sup>1</sup></a>	•	•	
<a href="#">Remote Plus app<sup>2</sup></a>	•		
<a href="#">uControl app</a>	•	•	
<a href="#">TV Connector<sup>2</sup></a>	•		
<a href="#">Remote control 2</a>	•	•	
<a href="#">uStream</a>	•	•	
<a href="#">uDirect 3</a>	•	•	
<a href="#">uTV 3<sup>3</sup></a>	•	•	•
<a href="#">uMic<sup>3</sup></a>	•	•	•
<a href="#">IIC remote (micro CIC remote)</a>	•	•	•
<a href="#">Remote control</a>			•
<a href="#">Smart Control remote</a>			•
<a href="#">uDirect 2</a>			•
<a href="#">uTV 2<sup>3</sup></a>			•
<b>For hearing healthcare professionals:</b>			
<a href="#">iCube II programming</a>	•	•	•
<a href="#">iCube programming</a>			•

<sup>1</sup>Product-specific rechargeable kit offered for Moxi All R, Moxi Fit R and Stride M R hearing instruments

<sup>2</sup>Compatible with Moxi All, 2.4 GHz hearing instruments

<sup>3</sup>uStream, uDirect 3 or uDirect 2 required

uControl app, uStream, uDirect 3, uTV 3, uMic, and IIC remote (micro CIC remote) are not compatible with the Moxi All hearing instruments  
Accessory compatibility may also be limited by product style and technology level.

Tempus platform:		North platform:		Era platform:	
T Max SP	T Stride P Dura	Stride M	Moxi Dura	Moxi <sup>2</sup>	Quantum <sup>2</sup> wireless ITEs
T Max UP	T Moxi Now	Stride P	Stride wireless ITEs	Moxi <sup>2</sup> Kiss	Quantum <sup>2</sup> non-wireless ITEs
T Moxi All	T Moxi Kiss	Stride P Dura	Stride non-wireless ITEs	Moxi <sup>2</sup> Dura	Max SPm
T Insera	T Moxi Fit	Moxi Now	Stride IIC	Quantum <sup>2</sup> S	Max SP
T Stride M	T Moxi Dura	Moxi Kiss		Quantum <sup>2</sup> HP	Shine Rev BTEs
T Stride P		Moxi Fit		Quantum <sup>2</sup> micro CIC	Shine Rev ITEs

## Accessory pairing guide

Product	uStream	uDirect 3	uDirect 2
<a href="#">Remote Plus app</a>			
<a href="#">uControl app</a>	Remote control	Remote control and My Ratings	No
<a href="#">uTV 3</a>	Streaming & charging	Streaming & charging	Streaming & charging
<a href="#">uMic</a>	Streaming	Streaming	Streaming
<a href="#">uTV 2</a>	Streaming	Streaming & charging	Streaming & charging