

At a glance

SoundCore		T Pro	T 800	T 700	T 600	T 500
SoundNav Automatic program identifies and classifies sounds across environments, with a focus on conversations.	Music	•	•			
	Noise	•	•	•		
	Conversation in a crowd	•	•	•		
	Conversation in a small group	•	•	•		
	Conversation in quiet	•	•	•	•	
	Conversation in noise	•	•	•	•	
	Quiet	•	•	•	•	
	Total environments	7	7	6	3	AutoMic
SpeechPro Provides optimal speech understanding and localization, even in the most challenging environments.	SpeechPro	•	SpeechZone™ 2	SpeechZone		
	Speech Locator	•	•	•		
	Speech Focus	•				
	Dynamic Spatial Awareness	•				
Sound Conductor Dynamically balances features to provide speech understanding, comfort and natural sound.	Speech enhancement	•	•	•	•	•
	Noise reduction	•	•	•	•	•
	Adaptive directionality	Multiband	Multiband	Multiband	Multiband	•
Spatial Awareness Makes it easier to localize sound providing a more natural, realistic listening experience.	Spatial Awareness	Dynamic	Personalized	•		
	Pinna Effect	•	•	•	•	•
Sound Stabilization						
	AntiShock 2	•	•	•	•	•
	Wind control	•	•	•	•	•
	Feedback manager	•	•	•	•	•
	Natural Sound Balance	•	•	•	•	•
Experience innovations						
	Log It All	•	•	•	•	•
	Patient Ratings	•	•	•	•	•
	Data logging	•	•	•	•	•
Flex						
	Flex:trial	•	•	•	•	•
	Flex:upgrade		•	•	•	•
Convenience						
	Rechargeable	•	•	•	•	•
	DuoLink	•	•	•	•	•
	Telecoil and Easy-t	•	•	•	•	•
	Binaural Phone	•	•	•	•	
Fitting						
	Automatic Adaptation Manager	•	•	•	•	•
	MyMusic	Binaural	Binaural	•	•	•
	Frequency compression	•	•	•	•	•
	Tinnitus masker	•	•	•	•	•
	IntelliVent	•	•	•	•	•
	Streaming programs	•	•	•	•	•
	Manual programs	•	•	•	•	•
	NAL-NL2/NL1 and DSLv5	•	•	•	•	•
	Fitting channels	20	20	16	10	6



Products shown at actual size
 Compatibility may be limited by product style