


# Stride at a glance

		Pro	800	700	600	500
Channels		20	20	16	10	6
<b>Signature features</b>	<b>Glossary</b>					
SpeechZone 2	In noisy situations, adjusts to an optimal microphone strategy to focus on speech from any direction	SpeechZone 2	SpeechZone			
Binaural spatial processing	Hearing instruments work together to determine the direction of speech	•	•			
SoundNav	Automatic program that identifies and smoothly transitions across environments, with a special focus on conversations	7 environments	6 environments	5 environments	2 environments	AutoMic
Sound Conductor	Dynamically balances the appropriate features providing speech understanding and comfort while maintaining natural sound quality	•	•	•	•	•
Log It All	Captures time spent across all seven SoundNav environments to support technology recommendations	•	•	•	•	•
MyMusic	Specialized fitter-controlled settings to enhance music enjoyment	Automatic binaural	Automatic binaural	•	•	•
Binaural Phone	Streams conversation from phone ear to opposite ear, no accessory required	•	•	•	•	
Automatic Adaptation Manager	Intelligently and gently transitions patients to the amplification required for long-term optimal speech understanding from first fit acceptance	•	•	•	•	•
<b>Features</b>						
Adaptive directional	Directional strategy that tracks and suppresses multiple background noise sources	Multiband	Multiband	Multiband	Multiband	•
Frequency compression	Compresses and shifts inaudible high-frequency sounds into a more audible lower frequency region	•	•	•	•	•
AntiShock	Eliminates discomfort from sudden impulse noises	•	•	•	•	•

## In all technology levels:

- Pinna Effect
- Natural Sound Balance
- Data logging
- Feedback manager
- Wind control
- Tinnitus masker
- Manual programs
- Streaming programs
- DuoLink
- Plasma coating
- Telecoil
- WDRC processing with NAL-NL1, NAL-NL2 and DSL v5
- Linear processing with DSL v5
- Easy-t



Accessories give you another tool to solve patient challenges. They are available at all technology levels.

